


Bioluminescence II – Serial Dilution with Luciferase *Protocol*


Reaction Set-up

1. Prepare a serial dilution of luciferase in your individual strip tubes. A serial dilution involves diluting a compound many times by taking a known concentration of a substance and moving a standard volume from one tube to the next in order to further dilute it.


- A. Pipette 55 μL from the tube labeled “luciferase” into the first tube (A) of your strip tube.



- B. Pipette 50 μL from the tube labeled “BUFFER” into tubes B - G of your strip tube.



- C. Pipette 5 μL of the luciferase solution from tube A and transfer it to tube B, pipetting up and down several times to mix.
- D. Repeat step C, but instead transfer 5 μL from tube B and transfer it to tube C, carefully pipetting up and down several times to mix. Continue transferring across the tubes (C→D, D→E, E→F) using the same technique.
- E. After you mix 5 μL into tube F, discard 5 μL from tube F into the waste container so that each tube now contains 50 μL liquid. **Do not transfer anything into tube G or H.**



2. Discuss with your group such that each of you will claim column 1, 2, 3, or 4 of your 96-well plate. Pick your columns in order and do not skip a column if your group is less than 4.
 - a. Write which column is yours here: _____

3. Using a multichannel pipette, move 40 μ l of each of your dilutions from your strip tube into your chosen column, rows A through G. Align your transfer such that tube A is transferred to row A and tube G is transferred to row G. Make sure each member of your group has added their dilution to the plate.

4. Using a multichannel pipette, move 40 μ l of Bright Glo reagent into Rows A \rightarrow H, Columns 1 \rightarrow 4.
 - a. Pour solution labeled "Bright Glo" into a plastic trough. Label your trough.
 - b. Use a multichannel pipette to add Bright Glo to rows A \rightarrow H of columns 1 \rightarrow 4.

	1	2	3	4	5	6	7	8	9	10	11	12
Row A	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row B	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row C	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row D	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row E	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row F	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row G	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								
Row H	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo	40 μ l Bright Glo								

View Results

After Bright Glo reagent is added, read your plate on a plate reader. Observe your data and write down any trends in amount of light output.

Questions

- A. Why did you perform your experiment in replicate (doing the same thing multiple times)?
- B. As you move from the top of the plate (Row A) to the bottom (Row H), does the amount of luciferase increase or decrease? Does it change by a similar amount or is it very different each time?
- C. As you move from the top of the plate (Row A) to the bottom (Row H), does the amount of light increase or decrease? Does it change by a similar amount or is it very different each time?
- D. When the data you generated was graphed so that the amount of luciferase in each well is on the x-axis and the corresponding amount of light produced on the y-axis, what was the shape of the data points?
- E. We already knew the concentrations of luciferase in each well after the serial dilution. Explain the use in creating a graph like the one you have generated.